

STX7-syntaxin 7 Gene

Reference: PTXBC011975

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	syntaxin-7; syntaxin 7		
For research use only.			

GENE DESCRIPTION	syntaxin 7
GENE ID	8417
ORIGIN SPECIES	Human
ACCESSIONS	BC011975
VECTOR	pDONR223
NCBI SYMBOL	STX7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttacctccaggagttgggtgacccgccagttggccagaggatcttcttaacatccagaagatcacacagtgttctgtggaatacaagaactctgaatcaacttgaacacctcaagattcacctgaattgaggcaacagttgcaacagaagcagcagtatactaccag
cttgccaaagaacagataagttacattaaagatttggatctctgccaccacccccagtgaaacagcgtcaaaaggaaaatacagaaggatcgttagtgccagagttcacaacatcactgacaaactccagaaggccagagcgtgctgagcagagagaagagttgtgtctcga
gtaagagccagttccagagtgctgagcagtttctgaggacagctcaaaagaaggaatctgtatcctgggaaagccaaactcaacctcaagtgacgtgaggtgaggaatgagaaattacagaggtgacctccgtcttattcatgagagagaatcttctatcaggcaactggaagctgatatta
tggatattaatgaaatatttaagatttgggaatgatgattcatgaacaaggagatgtaatatagatagcatagaagccaatgtggaaatgcagaggtgcacgttcagcaagcaaatcagcagctgtcaagggcagcagattatcagcgaatccagaaaacccctgtgcatcatcattctat
ccttgcattggagttgcgattatcagctcatcatatggggattgaaccactga