

CXCL1-chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha) Gene

Reference: PTXBC011976

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FSP; GRO1; GROa; MGSA; MGSA-a; NAP-3; SCYB1; growth-regulated alpha protein; C-X-C motif chemokine 1; GRO-alpha(1-73); GRO1 oncogene (melanoma growth stimulating activity, alpha); GRO1 oncogene (melanoma growth-stimulating activity); MGSA alpha; chemokine		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha)
GENE ID	2919
ORIGIN SPECIES	Human
ACCESSIONS	BC011976
VECTOR	pDONR223
NCBI SYMBOL	CXCL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccggcgctgctctccgccccagcaatccccgctcctgcgagtgccactgctgctctgctctgtagccgctggccggcgcgagcaggagcgtccgtggccactgaactgcgctgccagtcttcagaccctgcagggaattcaccccaagaacatccaaagtgtgaacgtgaagtcctcccgaccctcgcgccccaccgaagtcatagccacactcaagaatggcggaagcttgcctcaatcctgcatcccccatagttaagaaaatcatgaaaagatgctgaacagtgacaaatccaactga