

TAF12-TAF12 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 20kDa Gene

Reference: PTXBC011986

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TAF12 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 20kDa; TAF2J; TAFII20; transcription initiation factor TFIID subunit 12; TAFII-20/TAFII-15; TAFII20/TAFII15; TATA box binding protein (TBP)-associated factor, RNA polymerase II, J,		

For research use only.

GENE DESCRIPTION	TAF12 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 20kDa
GENE ID	6883
ORIGIN SPECIES	Human
ACCESSIONS	BC011986
VECTOR	pDONR223
NCBI SYMBOL	TAF12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaccagtttgccctcagccctaataacctccaattctcatcataaaacgggaaccagccagcaccctccacaaggtccatggccaatagtactgcagtggttaagataaccaggcactcctggggcaggaggtcgtcttagccctgaaaacaatcaggtattgaccaagaa
gaaattacaggacttagtaagagaagtgatcctaataagcagttggatgaagatgtggaggagatgctgctgcagattgctgatgatttatcgagagtggtgacagcagcctgtcagcttgcgcggcatcgcaagtctagcaccctggagggtgaaagatgtccagctgcatttagagc
gccagtggaacatgtggatcccaggatttggctctgaagaatccgaccctacaaaaagcttgcaccacagaagtcacaacagagaatggcattgatccggaaaacaaccaagaataa