

C5orf44-chromosome 5 open reading frame 44 Gene

Reference: PTXBC012006

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0533 protein C5orf44; C5orf44; trafficking protein particle complex subunit 13; Trs65-related; trafficking protein particle complex 13		

For research use only.

GENE DESCRIPTION	chromosome 5 open reading frame 44
GENE ID	80006
ORIGIN SPECIES	Human
ACCESSIONS	BC012006
VECTOR	pDONR223
NCBI SYMBOL	C5orf44
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagtgaatccccctaacaggagcacctgtggcgctaaagtgtgcggctgactaacctcttattaccaatataccagtaacatgtgaagagaaagacttacctggagatctcttaaccagctgatgagagatgatcctcaaccgttaagtgtgcagaagttaattgttggg
agaaatgctgactttaccacagaatttgggaatatatttgggagagacctttccagttatatcagcgttcataatgatagcaatcaagttgtaaaagacatattagttaaagctgatctcagacaaagttctcagcgtttaaatcttcagcctccaatgctgcagtggctgaacttaaccggatt
gtgtattgatgatgcatacatcatgaagtcaagaaattggaacacacatcttggtatgtgctgtgagtatacaactcaggctggagaaaaatgtattcagaaaaattctcaaatcaggttctcaaacattggatgtgaaaccaaattttacaatgcagagatgacctcagttctgtgac
tgatgaagtattctggaagccagattcagaataatgacaacctcacctatgtttatggagaaggtttcactggagccatctattatgtacaatgaacagaataaattcagtcagccaagctggagaatgtgtgtctacgtttgggtcaagagcatatttgaaccaatggatacagccagtac
ttatactgcctaaagccaaagaatga