

RPP14-ribonuclease P/MRP 14kDa subunit Gene

Reference: PTXBC012017

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HsHTD2; P14; ribonuclease P protein subunit p14; 3-hydroxyacyl-[acyl-carrier-protein] dehydratase; Hydroxyacyl-thioester dehydratase type 2, mitochondrial; ribonuclease P/MRP 14kDa subunit; ribonuclease P/MRP subunit p14		

For research use only.

GENE DESCRIPTION	ribonuclease P/MRP 14kDa subunit
GENE ID	11102
ORIGIN SPECIES	Human
ACCESSIONS	BC012017
VECTOR	pDONR223
NCBI SYMBOL	RPP14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgccctgctgccacatgaaagagtagttacaaaaaccctccgagtagccactacatgaaagtctgctagaattcaagattgtggagttggactgaatgctgcacagttcaaacagctgcttattcggctgtgaaggacctgttggggaggttgatgccgccttaccttggaca
tcctaacctatgaagagaagacctgtcagccatcttgagaatagtagcagtggtctgtgcaaatgtggagctcttgaccctgttagatcctataaaggcaaaaatgtgcttccgggtgattcaggtttctccatttctctgcattatctgtaatataggggaactagtattggattga