

RNF125-ring finger protein 125 Gene

Reference: PTXBC012021

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF125; TNORS; TRAC-1; TRAC1; T-cell RING activation protein 1; T-cell ring protein identified in activation screen; ring finger protein 125, E3 ubiquitin protein ligase; ring finger protein 125		

For research use only.

GENE DESCRIPTION	ring finger protein 125
GENE ID	54941
ORIGIN SPECIES	Human
ACCESSIONS	BC012021
VECTOR	pDONR223
NCBI SYMBOL	RNF125
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtcgtgctgagcaccacagcggcaaatcggcgccgctctgccaccgcggcgccctggagcgcaggaggacccggagttgccgtcacgtccttgcactgcgcgtgtgccttgagggtttacaccagcctgtccggaccgctgtggccacgtattctgccgttctgtattgctaccagctctaaagaacaacaagtgagactgtccttattgccgggcatactctcttcagaaggagttccagcaactgatgtagccaaaagaatgaaatcagagtataagaactgcgctgagtgtagacatagttctttacctcagtgaaatgagggcacatattcggactgtcagaagtacatagataaglatggaccactacaagaacttgaggagacagcagcaaggtgtgtatgtcccttttgcagaggggaactgtatgaagacagcttgcctggatcattgtattactatcacagatcggacggaggcctgtgttctgtccacttgcggtttaataccgatgagaatccaagcagcttcagtgccagtttaataagacatctgcaagttagtcacactttgtttatgatgattcatagatttaataataatgaggaagctcttatccgaagagcttagaccggtcacttcttgaatatgtgaatcactcgaacacagcataa