



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CENPH-centromere protein H Gene

Reference: PTXBC012024

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	centromere protein H; CENP-H; interphase centromere complex protein 35; kinetochore protein CENP-H		

For research use only.

GENE DESCRIPTION	centromere protein H
GENE ID	64946
ORIGIN SPECIES	Human
ACCESSIONS	BC012024
VECTOR	pDONR223
NCBI SYMBOL	CENPH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtagcagcagccccaatgatcaagacgcgcagcccgcgactcggaggggaagccggcgagccggcgccaccagcaggtgcggcgccccaagcggcggtgcagcagaggaccgcgatgccctgctctcaggctgagagcagacagaaaaacaactctgaatat
aaatcaatggttgatcgaagtgaagaaaaaacctcgaacaacaattatcgaagaaagcgaatcgaagctaaattgaagaccctggaaaatgaattgaagaggtaaaaagtgctgtttgagataaaaaagcttgcttagacaggaatgagactttcaactgcacttaaaaaaacctggagaaa
attagcagcagctctagtgtgctcattggataacatgaacacctattagagctcaataaataataatgaatcacagcaggaatctgggatttagagaaaacactgctgtatattagaagaagagagtgccaataaaacaagcttcgaagaagtgaagcttttagaataacagactgaaagaac
aaacagaaagtattggttcagacgatggaaaaactcagagagcagataaaatgatcatcagacaaaacctacagcttcgagatgaataaaaataactactctgtattcaacatgtgttcagaaccttatttggggagtaaaagtcattgggcagaggaatctgcccttaaggaaattgttctgcagcttgagaa
gaatttgacatgatgtaa