

GOLGA7-golgi autoantigen, golgin subfamily a, 7 Gene

Reference: PTXBC012032

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GOLGA3AP1; GOLGA7A; HSPC041; golgin subfamily A member 7; Golgi complex-associated protein of 16kDa; golgi autoantigen, golgin subfamily a, 7; golgi complex-associated protein of 16 kDa; golgin A7		

For research use only.

GENE DESCRIPTION	golgi autoantigen, golgin subfamily a, 7
GENE ID	51125
ORIGIN SPECIES	Human
ACCESSIONS	BC012032
VECTOR	pDONR223
NCBI SYMBOL	GOLGA7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggccgcagcaggcgcgggtgtccggaagggtgttcattcagcagactacagcagtcgacacgctgccagttccagacacagtcctcggagctggagaaccggattgataggcagcagtttgaagaacagttcgaactctaataacctttatgcagaagcagagaagct
tggcgccagtcatactcgaagggtgtttgctgttttaacgcatataccatctctatgcatggaaactcattatgagaaggttctgaagaaggtctcacaatacattcaagagcagaatgagaagatctatgctccacaaggcctcctctgacagaccctattgagcgaggactcgaggt
tattgaaattaccattatgaagacagaggcatgagcagtggaagataa