

HES2-hairy and enhancer of split 2 (Drosophila) Gene

Reference: PTXBC012091

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	bHLHb40; transcription factor HES-2; class B basic helix-loop-helix protein 40; hairy and enhancer of split 2; hes family bHLH transcription factor 2		

For research use only.

GENE DESCRIPTION	hairy and enhancer of split 2 (Drosophila)
GENE ID	54626
ORIGIN SPECIES	Human
ACCESSIONS	BC012091
VECTOR	pDONR223
NCBI SYMBOL	HES2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggctgcctcgcgggcaggggagcggcgagctgcgcaagcctgaagccgctgctggagaagcggcgcgcgcatcaaccagagcctgagccagcttaaggggctcatcctgccgctgctgggcccggaggatgcttctggtggcacacctggctccctc
catgctcagaactgcttctactctacatccaggctcctgagcagccccagcttaa