



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

NKIRAS1-NFKB inhibitor interacting Ras-like 1 Gene

Reference: PTXBC012145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KBRAS1; kappaB-Ras1; NF-kappa-B inhibitor-interacting Ras-like protein 1; I-kappa-B-interacting Ras-like protein 1; NFKB inhibitor interacting Ras-like protein 1; kappa B-Ras protein 1; kappa B-ras 1; NFKB inhibitor interacting Ras like 1		

For research use only.

GENE DESCRIPTION	NFKB inhibitor interacting Ras-like 1
GENE ID	28512
ORIGIN SPECIES	Human
ACCESSIONS	BC012145
VECTOR	pDONR223
NCBI SYMBOL	NKIRAS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr