

CD27-CD27 molecule Gene

Reference: PTXBC012160

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD27 molecule; T-cell activation antigen CD27; CD27 antigen; S152; S152. LPFS2; T14; TNFRSF7; Tp55; CD27L receptor; T cell activation antigen S152; tumor necrosis factor receptor superfamily, member 7		

For research use only.

GENE DESCRIPTION	CD27 molecule
GENE ID	939
ORIGIN SPECIES	Human
ACCESSIONS	BC012160
VECTOR	pDONR223
NCBI SYMBOL	CD27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcacggccacatccctgggtgctgtgcgttctggggaccctgggtgggctctcagctactccagcccccaagagctgccagagaggcactactgggtcagggaaagctgtgtgccagatgtgtgagccaggaaacattctcgtgaaggactgtgaccagcatagaaaggctg
ctcagtgatccttgataccgggggtctctctctcctgaccaccacacccggccccactgtgagagctgtcggcactgttaactctgttcttctgttcgcaactgacacatcactgccaatgctgagtgctgtcgaatggctggcagtgacaggacaaggagtgacccagtggtg
atctcttccaaaccttgcgtgaccgctcggtcgtctcagggcctgagccacacccctcagccacccacttaccttatgtcagtgagatgctggaggccaggacagctgggcacatgcagactctggctgacttcaggcagctgctgccccgactctcttaccactggccacccaa
agatccctgtgcagctcggatatttattcgcacctgtgatcttcttggaatgttctgttttcaacctggccggggccctgttctccatcaacgaagaaatatagatcaaaagagaaagtcctgtggagcctgcagagcctgtgttacagctgccccaggaggaggaggcag
caccatcccatcaggaggattaccgaaaaccggagcctgctgtccccctga