



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

THAP11-THAP domain containing 11 Gene

Reference: PTXBC012182

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CTG-B43a; CTG-B45d; HRIHFB2206; RONIN; THAP domain-containing protein 11; THAP domain containing 11		

For research use only.

GENE DESCRIPTION	THAP domain containing 11
GENE ID	57215
ORIGIN SPECIES	Human
ACCESSIONS	BC012182
VECTOR	pDONR223
NCBI SYMBOL	THAP11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr