

PLAC8-placenta-specific 8 Gene

Reference: PTXBC012205

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15; DGIC; PNAS-144; onzin; placenta-specific gene 8 protein; down-regulated in gastrointestinal cancer protein; placenta specific 8		

For research use only.

GENE DESCRIPTION	placenta-specific 8
GENE ID	51316
ORIGIN SPECIES	Human
ACCESSIONS	BC012205
VECTOR	pDONR223
NCBI SYMBOL	PLAC8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaagctcaggcgccggtgtgtgtgacccaacctggagtcgggtccggtccgccccagaaacccaactggcagacagcatgtgtgactgtttcagcgactcggagtcgtctctgtgacacattttttcccggtcccttggtgcaagttgcagctgatataatgaatgctgtctgtgtggaacaggcgtcgcaatgaggactctctacaggaccgatatggcatccttgatctatttgtgactatattggcaactctttgctgtctcattgtactctttgccaatcaagccatctgattacagaaggagaccatgcgtactttctaa