

ACYP2-acylphosphatase 2, muscle type Gene

Reference: PTXBC012290

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACYM; ACYP; acylphosphatase-2; acylphosphatase 2, muscle type; acylphosphatase, muscle type isozyme; acylphosphate phosphohydrolase 2; testicular tissue protein Li 11; acylphosphatase 2		

For research use only.

GENE DESCRIPTION	acylphosphatase 2, muscle type
GENE ID	98
ORIGIN SPECIES	Human
ACCESSIONS	BC012290
VECTOR	pDONR223
NCBI SYMBOL	ACYP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctaccgccagtcactcaaatccgtggactacgaggtgttcgaagagtgacgggtgtttgcttcagaatgtatacagaagatgaagctaggaaaataggagtggttggtggtgaagaataccagcaaaaggcaccgtgacaggccaagtgcaggggccagaagacaaagtcaa
ttccatgaagtcctggctgagcaaggttggaagccctagttctgcattgaccgcacaaacttttctaataaaaaaacatctctaagcttgaatactctaatttttagtattagatactaa