

RPL39L-ribosomal protein L39-like Gene

Reference: PTXBC012328

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RPL39L1; 60S ribosomal protein L39-like; 60S ribosomal protein L39-2; L39-2; ribosomal protein L39-like 1; ribosomal protein L39-like protein; ribosomal protein L39 like		

For research use only.

GENE DESCRIPTION	ribosomal protein L39-like
GENE ID	116832
ORIGIN SPECIES	Human
ACCESSIONS	BC012328
VECTOR	pDONR223
NCBI SYMBOL	RPL39L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttcacaaagacttcaccattagcgattcctggccaagaacaaaagcaaaatcgccatccccagtggttcagatgaaacctggttagtaaaatcaggtacaactccaaaaggagcattggagaagaaccaagctgggtctataa