



Reference: PTXBC012349

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBX1; FBXO1; cyclin-F; F-box only protein 1; G2/mitotic-specific cyclin-F; cyclin F		
For research use only.			

GENE DESCRIPTION	cyclin F
GENE ID	899
ORIGIN SPECIES	Human
ACCESSIONS	BC012349
VECTOR	pDONR223
NCBI SYMBOL	CCNF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr