



Reference: PTXBC012359

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLPTM1, transmembrane protein; cleft lip and palate transmembrane protein 1; cleft lip and palate associated transmembrane protein 1		

For research use only.

GENE DESCRIPTION	cleft lip and palate associated transmembrane protein 1
GENE ID	1209
ORIGIN SPECIES	Human
ACCESSIONS	BC012359
VECTOR	pDONR223
NCBI SYMBOL	CLPTM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr