

C2orf56-chromosome 2 open reading frame 56 Gene

Reference: PTXBC012374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C2orf56; MidA; PRO1853; protein arginine methyltransferase NDUFAF7, mitochondrial; NADH dehydrogenase (ubiquinone) complex I, assembly factor 7; NADH dehydrogenase [ubiquinone] complex I, assembly factor 7; mitochondrial dysfunction protein A homolog; pro		

For research use only.

GENE DESCRIPTION	chromosome 2 open reading frame 56
GENE ID	55471
ORIGIN SPECIES	Human
ACCESSIONS	BC012374
VECTOR	pDONR223
NCBI SYMBOL	C2orf56
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcacagggaagtagcctctgggcccaataaaacacacacatttttaaaataggggtattgatgtccggctgaaggtcttttagataaatcaaatgagccatcagtgaggcagcagttactcaaggatatgatgttaatgaatccaagaagatgggagagagatttaacttttgccttgcctacatcagagactcaaggtggaagatcagagggaatgcagtcagtcacaaaccccttgcacccgtttagctgggttagtgaacttgcctggcagtgga