

PDCD6-programmed cell death 6 Gene

Reference: PTXBC012384

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALG-2; ALG2; PEF1B; programmed cell death protein 6; apoptosis-linked gene 2 protein; programmed cell death 6		

For research use only.

GENE DESCRIPTION	programmed cell death 6
GENE ID	10016
ORIGIN SPECIES	Human
ACCESSIONS	BC012384
VECTOR	pDONR223
NCBI SYMBOL	PDCD6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgctactctaccgccccggcctggggccggcctgggcctgtgcaggcgcgcgtgccggaccagagcttctgtgaacgtttccagaggggtcgataaagacaggagtggagtgtatcagacaccgagcttcagcaagctctctccaacggcacgtggactccctt
taatccagtgactgtcagggtcgtatcattatccatgttgaccgtgagaacaaggccggcgtgaacttcagcaggttcacgggtgtgtggaagtacatcacggactggcagaacgtctccgcacgtacgacgggacaactccgggatgatcgaataagaacgagctgaagcaggccctctc
agggttcgggtaccggtctctgaccagttccacgacatcctcattcgaaggtttgacaggcaggggacggggcagattgccttcgacgacttcacccagggtgcacgtcctgcagaggttgacggatatattcagacgttacgacacggatcaggacgggtggttcaggtgtcgtacg
aacagtacgtgtccatggtcttcagtatcgtatga