

FBXO17-F-box protein 17 Gene

Reference: PTXBC012385

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBG4; FBX26; FBXO26; Fbx17; F-box only protein 17; F-box only protein 26; F-box protein FBG4; F-box protein 17		

For research use only.

GENE DESCRIPTION	F-box protein 17
GENE ID	115290
ORIGIN SPECIES	Human
ACCESSIONS	BC012385
VECTOR	pDONR223
NCBI SYMBOL	FBXO17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggcccccgtatcgcgcgacggctgccggcggaccgtccctggccctggacgcgctccccggagctgctggtgagccacgtccgccacgtcttggtcacgcgatgccgccagtgccgcctggcgcgacataggacggggccactgtgtggctgctgacgtggcccgacgcgacgcgcggagggccgcgacgtctacgagtggtcaacgtgctgccccagcaacgaagacaaggaggtcccgctgtgcccctggcgcgctactgtctgcccgcgcccctggccgcaatctcatcttcaactcctgcggagagcaggcttcagaggctgggaggtggagcatggcgggaacggctggccatagaaaagaacctaacacgggtgctgggctccttcgcagacctgctctgtgacctcttgcgaatgggtgctccaaggagggcagcttggacctggtgatggaaggggtgtggcaggagctgctggacagcggccagattgagatctgtgtggctgactggtggggcgctcgagagaactcgggctgctaccagctccgggtccgcttctggtatgtatgaaaaggagtggtcaagttctcagctcacctgacccggctcctcagtgactgagaggggtgccgacaggtctccacgtcttccaaactttgcaaggcatccgctacgtatctttgagcagtacgggagagacgtgagttcctgggtggggcactacggcccttgtgaccactccagtgtaggggtcaggatccgtctgtcctag