



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

WRB-tryptophan rich basic protein Gene

Reference: PTXBC012415

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tail-anchored protein insertion receptor WRB; CHD5; GET1; congenital heart disease 5 protein; tryptophan rich basic protein		

For research use only.

GENE DESCRIPTION	tryptophan rich basic protein
GENE ID	7485
ORIGIN SPECIES	Human
ACCESSIONS	BC012415
VECTOR	pDONR223
NCBI SYMBOL	WRB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctcaccgcccggccgaccattggcgtggttctgctctcagcttctggtttgagtcgaatgcttctaggtacctctccgcctctctctcatctctcatgtccaggtgctctgcagaaaggacgaggagcaggtgtcacagatgagagcggagatcaggacatgaagcaggagctctccaatgtcaacatgattggagacgtatttggccagatgtccgagctggcgaagaagaatcacacagatgacgtagatcaactcaaaaacatgtgaaagctgcagacatcaattagcaggaataaaatgggtgtgaaatgtgtctctcattgtattgtcagatgtcagagctgcctgatgatctcaactatgtgaattatttcttccctctgtctctgtctctgcagataaaatgtcaaacccctcagaccctgtgagcttctctactagatgacacgtgtgtgaaatcttctgttattgttgcattctcaaacaaatctgcattatgtctctcaatcactgac