



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD40-CD40 molecule, TNF receptor superfamily member 5 Gene

Reference: PTXBC012419

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD40 molecule; CD40 molecule, TNF receptor superfamily member 5; B cell surface antigen CD40; Bp50; CDW40; TNFRSF5; p50; tumor necrosis factor receptor superfamily member 5; B cell-associated molecule; CD40L receptor		

For research use only.

GENE DESCRIPTION	CD40 molecule, TNF receptor superfamily member 5
GENE ID	958
ORIGIN SPECIES	Human
ACCESSIONS	BC012419
VECTOR	pDONR223
NCBI SYMBOL	CD40
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]