

YIPF6-Yip1 domain family, member 6 Gene

Reference: PTXBC012469

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein YIPF6; FinGER6; YIP1 family member 6; Yip1 domain family member 6		
For research use only.			

GENE DESCRIPTION	Yip1 domain family, member 6
GENE ID	286451
ORIGIN SPECIES	Human
ACCESSIONS	BC012469
VECTOR	pDONR223
NCBI SYMBOL	YIPF6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaagcggaggagtcctcaggagacccgggacagcatgccagccctgtttgcaggccttcagatatcatctcacaagacatcccgtagaaggagaaatcaccattcctatgagatctcgcatcgggagtttgacagctccacattaaatgaatctgttcgaatac
catcatgcgtgatctaaaagctgttgggaaaaattcatgcatgtttgtacccaaggaaaagtaatactctttgagagattgggatttggggccctttgatcccttggtagacactcgcattaatgctgcaaaagagactctgagatagtgaaaaagatggaggccccaattgcagagggtg
tttgcattgtctgttgggtgagttaccatcacctcaactcaaaactcttggagggaacafatctttttcagagccctctgtgtgctgggttactgtatacttcccttgacagtagcaatgctgattggcggctggtacttttggctgatccaggacctgtaaaactcatggttcggcttttgggt
gattgtgatgttgcctggtctatagttgcctccacagctctcctgtgatagccagcctcaaacccagagccctagctgtttatcctgttttctgtttactttgtcatcagttggatgattcacccttactcctcagtaa