

CLDN1-claudin 1 Gene

Reference: PTXBC012471

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLD1; ILVASC; SEMP1; senescence-associated epithelial membrane protein 1; claudin 1		

For research use only.

GENE DESCRIPTION	claudin 1
GENE ID	9076
ORIGIN SPECIES	Human
ACCESSIONS	BC012471
VECTOR	pDONR223
NCBI SYMBOL	CLDN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaacgcgggctgcagctgttgggttcattctgccttctgggatggatcggccatcagcactgccctgccccagtgaggattactctatgccggcgacaacatcgtgaccgccaggccatgtacgagggctgtggatgtcctgcgtgtcgcagagcaccggg
cagatccagtgcgaagtcttggactccttgcctgaatctgagcagcacttgcgaagcaaccgtgcttgatgggttggtggcatctcctgggagtgatagcaatcttggccaccgttgcatgaagtgtatgaagtcttggagacgatgaggtgcagaagatgaggatggctgtcattg
ggggcgcatattcttcttgcaggtctggctattttagtgcacagcatggtatggcaatagaatcgttcaagaattctatgacctatgacccagtcgaatggcagggtacgaatttggcaggctcttctactggctgggctgcttctctgcttcttgggaggtgcctactttgctgttc
ctgtccccgaaaaaacctcttaccacaacacagccctatccaaaacctgcacctccagcgggaaagactacgtgtga