

UBD-ubiquitin D Gene

Reference: PTXBC012472

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UBD-3; FAT10; GABBR1; ubiquitin D; diubiquitin; ubiquitin-like protein FAT10		
For research use only.			

GENE DESCRIPTION	ubiquitin D
GENE ID	10537
ORIGIN SPECIES	Human
ACCESSIONS	BC012472
VECTOR	pDONR223
NCBI SYMBOL	UBD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctccaatgcttctgtgcatgtccgtccgaggaatgggatttaacacccatgatgacagcgtgaaaaaaatcaagaacatgtccgtctaaagccaaggttctgtgaggaccaggtctttgtgggtccaagatctaaagccacggagaag
cctctcatcttatggcattgacaaagagaagaccatccactacccctgaaagtggggaagccagtgatgaggagctgccctgtttctgtgagtgaggtgatgaggcaagagagcactcctccaggtgcgaaggtccagctcagtgccacaagtgaagcaatgatcgagactaag
acgggtataatccctgagaccagattgtgactgcaatggaaagagactggaagatgggaagatgatggcagattacggcatcagaagggaacttactcttctgcatgttattgattggagggtga