

KIAA1147-KIAA1147 Gene

Reference: PTXBC012493

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LCHN; PRO2561; protein LCHN		

For research use only.

GENE DESCRIPTION	KIAA1147
GENE ID	57189
ORIGIN SPECIES	Human
ACCESSIONS	BC012493
VECTOR	pDONR223
NCBI SYMBOL	KIAA1147
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttcttgaattgattttgtctgtgtacaaataacatctatgaatcatcatctgtatatctgaaccaagttgtgtgcagagaggttgatataactatattacagaaatgattttacaattaa