

POP5-processing of precursor 5, ribonuclease P/MRP subunit (*S. cerevisiae*) Gene

Reference: PTXBC012505

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POP5 homolog, ribonuclease P/MRP subunit; ribonuclease P/MRP protein subunit POP5; HSPC004; RPP2; RPP20; hPop5; processing of precursor 5, ribonuclease P/MRP subunit		

For research use only.

GENE DESCRIPTION	processing of precursor 5, ribonuclease P/MRP subunit (<i>S. cerevisiae</i>)
GENE ID	51367
ORIGIN SPECIES	Human
ACCESSIONS	BC012505
VECTOR	pDONR223
NCBI SYMBOL	POP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgcggttcaagcacaggtacgtgctctgcaactggtgtctgacgacccccgctgccgctaagcctcgatgaccgagttctgagcagcctctgacgggacacgacgccaggggtgcacggaactttcggcgagccgctgctccatggcttcggttcgatatctcaatgct
atactggaatagtgctacttcgatgcagaaaagaattctatcagcttgtgtgtcagctctcccttcacatacttgagagaacaaaggacaccgttacccatgcttttcaacacattacatgtgggaggtacaataagaacatgtcagaagttcctaattcagtaacacaggagacagctgttg
atctgtgtcagaactgcactgatgaaggagagcgggaagctatccagaagctgtgacaagaagctgcttactagaggaggaggagtcagggtgaggaggctgcagaagcaatggagtga