

IL11-interleukin 11 Gene

Reference: PTXBC012506

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGIF; interleukin-11; adipogenesis inhibitory factor; oprelvekin; interleukin 11		
For research use only.			

GENE DESCRIPTION	interleukin 11
GENE ID	3589
ORIGIN SPECIES	Human
ACCESSIONS	BC012506
VECTOR	pDONR223
NCBI SYMBOL	IL11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaactgtgttgccgctggctcctggtcgtgctgagcctgtggccagatacagctgtcggccctgggcccacacctggccccctcagtttcccagaccctgggcccagctggacagaccgtgctcctgacctctctctcctggcgacacgcggcagctggctgcacagctg
agggacaaattcccagctgacggggaccacaacctggattccctgcccacccctggccatgagtcgggggacactgggagctctacagctcccaggtgtgtgctgacaaggctgcgagcgggacactgtcctacctgcggcagctgagtggtgctgcggggcaggtggtcttccctga
agacctggagcccagctgggacacctgcaggcccgactggaccggctgctgcggcgtgagctcctgatgtccgcctggccctgccccagccacccccggaccgccggcgccccctggtggccccctctcagctggggggcatcaggggccgcctcgcacatc
tgggggggctgcacctgacacttgactggggcctgagggggactgctgctgctgaagactcggctgtga