

CXCL11-chemokine (C-X-C motif) ligand 11 Gene

Reference: PTXBC012532

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H174; I-TAC; IP-9; IP9; SCYB11; SCYB9B; b-R1; C-X-C motif chemokine 11; beta-R1; chemokine (C-X-C motif) ligand 11; interferon gamma-inducible protein 9; interferon-inducible T-cell alpha chemoattractant; small inducible cytokine B11; small inducible cyto		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) ligand 11
GENE ID	6373
ORIGIN SPECIES	Human
ACCESSIONS	BC012532
VECTOR	pDONR223
NCBI SYMBOL	CXCL11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgtgaaggcatggctatagccttgctgtgatattgtgtctacagttgttcaaggcttccccgtgtcaaaagaggacgctgtcttgcataggccctgggttaaaagcagtgaaagtggcagatattgagaaagcctccataatgtacccaagtaacaactgtgacaaaatagaagt
gattattaccctgaaagaaaataaaggacaacgatgcctaaaccccaaatcgaagcaagcaaggcttataatcaaaaaagttgaaagaagaatttttaa