

MED31-mediator complex subunit 31 Gene

Reference: PTXBC012539

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	3110004H13Rik; CGI-125; Soh1; mediator of RNA polymerase II transcription subunit 31; hSOH1; mediator complex subunit SOH1; mediator of RNA polymerase II transcription, subunit 31 homolog; mediator complex subunit 31		

For research use only.

GENE DESCRIPTION	mediator complex subunit 31
GENE ID	51003
ORIGIN SPECIES	Human
ACCESSIONS	BC012539
VECTOR	pDONR223
NCBI SYMBOL	MED31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgctgctgctatggagacagatgatgctggaatcgactcgggttcagttggagttggaattgtgcaatgttagccaacccaattacctaattttctgccaaagaggttactcaagacaaagctttgttaattatctaaatacttgccttactggaagaccagaatatgc
caagtatctaaagtaccctcagtggtttacacatgttagagctgctccaatgatgaacacttcgaaaggagctggtgaatgctcagtgctgcaatttattgatgaacagcagattctacattggcagcactattccggaagcggatgcgcttcagcaagccttggcagagcagcaacagcaa
aataacacatcgaggaaatga