

ORMDL2-ORM1-like 2 (*S. cerevisiae*) Gene

Reference: PTXBC012543

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC160; MST095; MSTP095; adoplín-2; ORM1-like protein 2; expressed in normal aorta; ORMDL sphingolipid biosynthesis regulator 2		

For research use only.

GENE DESCRIPTION	ORM1-like 2 (<i>S. cerevisiae</i>)
GENE ID	29095
ORIGIN SPECIES	Human
ACCESSIONS	BC012543
VECTOR	pDONR223
NCBI SYMBOL	ORMDL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgtgggggtggcacacagcgaagtaaaccccaacacccgagtgatgaatagccgagcatctggctggcctacatccttggtaggattgctgcatatgggttctactcagcatccctctcagcattcctgtgtctggaccctgaccaacgtcatcataacctggctacgtatgtcttcctcaccgtgaaaggacaccccttgagactcctgaccaaggaaagcccggtactgacacactgggagcaaatggactatgggtccagtttacctcttcccgaagttctcagcagctctcctattgtgctctatctcctggccagcttctataccaagtatgatgctgcgcacttctcatcaacacagcctcattgctaagtgtactgtgcccgaagttgcccagttccatggggttcgtgtcttggcatcaacaatactga