

PLUNC-palate, lung and nasal epithelium associated Gene

Reference: PTXBC012549

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein Plunc; PLUNC; LUNX; NASG; SPLUNC1; SPURT; bA49G10.5; BPI fold-containing family A member 1; ligand-binding protein RYA3; lung-specific protein X; nasopharyngeal carcinoma-related protein; palate lung and nasal epithelium clone protein; palate, lun		

For research use only.

GENE DESCRIPTION	palate, lung and nasal epithelium associated
GENE ID	51297
ORIGIN SPECIES	Human
ACCESSIONS	BC012549
VECTOR	pDONR223
NCBI SYMBOL	PLUNC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttcaaacctggggcctcattgtctctacgggctgttagccagaccatggccagtttgaggcctgccctggaccagacctgccctgaatgtgaatccagcctgccctgagtcacaggtctgcaggaagcttgacaaatgccctcagcaatggcctgctgtctg
ggggcctgttgggcatcttggaacacccctcctctgacatcctgaagcctggaggaggtactctgtggcctccttgggggactgcttggaaggtgacgtcagtgatcttgccctgaacacacattgacataaagtgactgacccccagctggaactggcctgtgcaga
ggcctgatggccaccgtctctatgtcaccatccctctggcataaagctcgaatgaatacggccctgtcgggtgcaagctgttgaggctggctggaagctggacatcactgcagaaatcttagctgtgagagataagcaggagagatccacctggtccttggctgacacccattccc
ctggagcctgcaaatctctgtctgactgagcctggccctcccatcaaggtctctggacagcctcacaggatctgaataaagctcctgagttggtcaggcgaacgtgtgccctctgtcaatgaggtctcagagcctggacatcacctggtgcatgacattgtaacatgct
gatccacggactacagttgtcatcaaggctaa