

MRPS10-mitochondrial ribosomal protein S10 Gene

Reference: PTXBC012560

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MRP-S10; PNAS-122; 28S ribosomal protein S10, mitochondrial; S10mt; mitochondrial 28S ribosomal protein S10; mitochondrial ribosomal protein S10		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S10
GENE ID	55173
ORIGIN SPECIES	Human
ACCESSIONS	BC012560
VECTOR	pDONR223
NCBI SYMBOL	MRPS10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggcggcagcgttcggtgctgtgtccggcgcctctggcagggattggggaattttctgtaaacacttctaaggcgaatacagccaaaatggcgttgccttcagtaaccaatgaagtggtacagtttcaaacctacacgttgatgttccaaaggattgacaaacctgt
ggtaacaatctctgatgaaccagacatattataaagcgcctctcgggtttgggtaagggtcacgataaggctgtattggacagttatgaatatttgcctgtcgtgctgtaaaagaacttgatctctattaaagtacatgaacctccaaggaataagagcgaattactcttccaatcagtgcat
atttacaagaagcacagaggttcagatgaaatgagaacactttacagatgttttagagttagaacatctaactggaagcacagcagatgtctacttggatattatcagcgaacttacctgaaggggttgccttggaagtaacaaagacacaattagaacagttaccagaacacatcaaggagc
caatctgggaacactatcagaagaaaaagaaagcaagtcataa