

CXCL13-chemokine (C-X-C motif) ligand 13 Gene

Reference: PTXBC012589

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANGIE; ANGIE2; BCA-1; BCA1; BLC; BLR1L; SCYB13; C-X-C motif chemokine 13; B-cell chemoattractant; B-cell-attracting chemokine 1; B-cell-homing chemokine (ligand for Burkitt's lymphoma receptor-1); B-lymphocyte chemoattractant; CXC chemokine BLC; b cell-at		

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) ligand 13
GENE ID	10563
ORIGIN SPECIES	Human
ACCESSIONS	BC012589
VECTOR	pDONR223
NCBI SYMBOL	CXCL13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagttcatctcgacatctctgcttctcatgctggtcagcagcctctctccagtcgaaggtgtctgaggtctattacacaagcttgaggtgtagatgtgtccaagagagctcagctttatccctagacgcttcattgatgaattcaaatcttgccccgtgggaatggtgtccaagaaaa
gaaatcatagctcgaagaagaacaagtcattgtgtgtgtgaccccaagctgaatggatacaagaatgatggagatttgagaaaaagaagttcttcaactctaccagttccagtggtttaagagaagaagattccctga