

TNNI1-troponin I type 1 (skeletal, slow) Gene

Reference: PTXBC012600

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SSTNI; TNN1; troponin I, slow skeletal muscle; troponin I type 1 (skeletal, slow); troponin I, slow-twitch isoform; troponin I1, slow skeletal type		

For research use only.

GENE DESCRIPTION	troponin I type 1 (skeletal, slow)
GENE ID	7135
ORIGIN SPECIES	Human
ACCESSIONS	BC012600
VECTOR	pDONR223
NCBI SYMBOL	TNNI1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggaagtcgagagaaacccaagatcactgcctccgcaaaactcttgctgaagagcctgatgctgccaaggccaaggaatgctgggagcaggagcagaggagcgcgaggctgagaaggtgcgctacctggcagagcgcacatccccacgtgcagaccgtggcctgtccc
tcagtgcctgcaggacctgtgcccggagctgcacccaaggtggaggtggtgatgagagcgatacgcacattgaggccaatgctccacaacaccaggagattaaggacctgaagctgaaggtgatggacctccgtgggaagtcaagcggcccccgtgcgtcaggtccgt
gtctcggctgacgcatgctccggccctgctgggtccaagcacaaggtgtccatggatctgcgggccaacctcaagtctgtgaagaaggaagacacagagaaggagcggcctgtggaggtgggtgactggaggagaacgtggaggccatgtctggcatggaaggccggaaga
agatgtttgatgccccaagtctccgacctcacaatag