

GOSR1-golgi SNAP receptor complex member 1 Gene

Reference: PTXBC012620

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GOLIM2; GOS-28; GOS28; GOS28/P28; GS28; P28; Golgi SNAP receptor complex member 1; 28 kDa Golgi SNARE protein; 28 kDa cis-Golgi SNARE p28; Golgi SNARE 28 kDa; cis-golgi SNARE; golgi integral membrane protein 2		

For research use only.

GENE DESCRIPTION	golgi SNAP receptor complex member 1
GENE ID	9527
ORIGIN SPECIES	Human
ACCESSIONS	BC012620
VECTOR	pDONR223
NCBI SYMBOL	GOSR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcaggaccagcagttactgggaagatctcaggaaacaggtcgacagctggaaaaatgaacttgacctgaaactagtttccttcagcaaaactgtacaagttacagtcatagcagtacccgagatggaagacgcgacaggtatgtctgatacaacaccccttttaaatggatc
aagccaagacagaatgtttgagacaatggcgattgagattgaacaacttttggcaaggcttacaggggttaattgataaaatggcagaatataccaacagtgccaggtgtccctccttgatgcagccctgatgcatacattacagcgccatagagacataatgcaggtaatatattgggctcga
gatgttttattatcacaggagtttgggttttttttaacccctgattggatattgacatatatttaaaaggacaaagagaaaaatctgagaaaaataatgccatgctaaaggagactgttttacttttttgggggtcaccattataaaattctga