



Reference: PTXBC012745

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RRP12-like protein; KIAA0690; ribosomal RNA processing 12 homolog		
For research use only.			

GENE DESCRIPTION	ribosomal RNA processing 12 homolog (S. cerevisiae)
GENE ID	23223
ORIGIN SPECIES	Human
ACCESSIONS	BC012745
VECTOR	pDONR223
NCBI SYMBOL	RRP12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr