



Reference: PTXBC012806

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SMVT; sodium-dependent multivitamin transporter; Na(+)-dependent multivitamin transporter; Na+-dependent multivitamin transporter; solute carrier family 5 (sodium-dependent vitamin transporter), member 6; solute carrier family 5 (sodium/multivitamin and i		

For research use only.

GENE DESCRIPTION	solute carrier family 5 (sodium-dependent vitamin transporter), member 6
GENE ID	8884
ORIGIN SPECIES	Human
ACCESSIONS	BC012806
VECTOR	pDONR223
NCBI SYMBOL	SLC5A6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr