

VAR5-valyl-tRNA synthetase Gene

Reference: PTXBC012808

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G7A; VARS1; VARS2; valine--tRNA ligase; protein G7a; valRS; valine tRNA ligase 1, cytoplasmic; valyl-tRNA synthetase 2; valyl-tRNA synthetase		

For research use only.

GENE DESCRIPTION	valyl-tRNA synthetase
GENE ID	7407
ORIGIN SPECIES	Human
ACCESSIONS	BC012808
VECTOR	pDONR223
NCBI SYMBOL	VAR5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr