

## HYAL3-hyaluronoglucosaminidase 3 Gene

Reference: PTXBC012892

|                |                                                                                              |                     |                |
|----------------|----------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                          | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                         | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | HYAL-3; LUCA-3; LUCA3; hyaluronidase-3; lung carcinoma protein 3; hyaluronoglucosaminidase 3 |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | hyaluronoglucosaminidase 3                                                                       |
| <b>GENE ID</b>          | 8372                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC012892                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | HYAL3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgaccacgcaactgggccagccctggtgctgggggtggccctgtgctggtgtggccagccctaccacaggtccctgaacgccccttctctgtgctggaatgtacctcagcacactgtgagggccctgttgggtgacactgccactcaatgctctgggcatcatagccaacc  
gtggccagcattttcacgggtcagaacatgacattttctacaagaaccaactggccctctatccctactttggaccaggggcacagctcacaatggggcatccccaggtttgccccttgaccgccacctggcactggctgctaccagatctaccacagcctgagacctggcttctg  
ggccagcagtgctggattggaggagtggtgtccactctgggctgggaactggggccgcccagcttatacaggcagcctcttgggcttggggcacagcaggtattccctgacctggacctcaggagcagctctacaaggcctatactggctttagcagggcggccctgcaactgatg  
gaggatacgtgctgggtggccaggcactacggcccatggactctggggcttctactaccacagcctgtggcaatggctggcatagtgatggcttcaactataccggccgctgcatgcagccacctggccgcaacactcaactgattggctctggggccctccagtgccctct  
tccccagcatctacctccccagggctgccacctgcccaccagggcctttgtccgacatgcctggaggaggccttccgtgtggccctgttggggcaccgacatccccctgctgctggtgctatgtccgctcacacaccggagatctgggaggttccgtgtccaggaatgacctgtg  
cagtccattggtgtgagtgacgactaggggcagccggcgtgtgtgctctgggggacctgagcctctccagctctgaggaggagtgctggcatctccatgactacctgtgtggacacctggggccctatgtgatcaatgtgaccagggcagcgtatggcctgcaatcaccagcgtgccca  
tggccacggggcgtgtgcccggcgagatccaggacagatggaagcctttctacacctgtggccagacggcagccttgagattggaagtccttcagctgccactgttactggggctgggctggccccacctgcccaggagccagcctgggctaaagaagcagtataa