

S100A3-S100 calcium binding protein A3 Gene

Reference: PTXBC012893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S100E; protein S100-A3; S100 calcium binding protein A3		
For research use only.			

GENE DESCRIPTION	S100 calcium binding protein A3
GENE ID	6274
ORIGIN SPECIES	Human
ACCESSIONS	BC012893
VECTOR	pDONR223
NCBI SYMBOL	S100A3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaggcctctggagcaggcggtagctccatcgtgtgcaccttcagggaatacgcagggcgctgtgggacaatacaagctctgccaggcggagctcaaggagctgctgcagaaggagctggccacctggaccccgactgagttcgggaatgtgactacaacaattcatga
gtgttctggacaccaacaaggactgcgaggtggacttftgtggagtatgtgcctcacttgcctctgtcttactgccacgagtacttcaaggactgcccctcagagccccctgtctccagtag