

SNAI1-snail homolog 1 (Drosophila) Gene

Reference: PTXBC012910

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein SNAI1; SLUGH2; SNA; SNAH; SNAIL; SNAIL1; dJ710H13.1; protein sna; protein snail homolog 1; snail 1 homolog; snail 1 zinc finger protein; snail family zinc finger 1; snail homolog 1; snail family transcriptional repressor 1		

For research use only.

GENE DESCRIPTION	snail homolog 1 (Drosophila)
GENE ID	6615
ORIGIN SPECIES	Human
ACCESSIONS	BC012910
VECTOR	pDONR223
NCBI SYMBOL	SNAI1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcgctcttctcgtcaggaagccctccgacccaatcggaagcctaactacagcgagctgcaggactctaatacagagttacctccagcagccctacgaccaggccacctgctggcagccatccacctccggagatcctcaacccaccgctcgtccaatgctcatct
gggactctgtctggcgccccaagccagccaatgctgggctctccttcggctccaggagagtgccagggtggcagagctgacctcctgtcagatgaggacagtgggaaaggctccagccccagccccaccctaccggctccttcctctcctcactcagctctcttccttg
gaggccgaggcctatgctgcttccaggcttggccaagtgcaccaagcagctggccagctctctgaggccaaggatctccaggtcgaaaggccttcaactgcaatactgcaacaaggaatacctcagcctgggtgcccctcaagatgcacatccgaagccacacgtgccctgcgt
ctgcggaaacctgcgggaaggcctctctagggcctgggtgctacaaggccatgtccggaccacactggcgagaagccctctctgtccccactgcagccgtgctctcgtgaccgtccaaacctgcggggccacctccagaccactcagatgtcaagaagtaccagtgccaggcgt
gtgctcggaccttctccgaatgctcctgctccacaagcaccaagagtcgggtgctcaggatgtccccgtga