



Reference: PTXBC012932

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H1VPS45; SCN5; VPS45A; VPS45B; VPS54A; VSP45; VSP45A; vacuolar protein sorting-associated protein 45; h1Vps45; leucocyte vacuolar protein sorting 45; vacuolar protein sorting 45A; vacuolar protein sorting 45 homolog		

For research use only.

GENE DESCRIPTION	vacuolar protein sorting 45 homolog (S. cerevisiae)
GENE ID	11311
ORIGIN SPECIES	Human
ACCESSIONS	BC012932
VECTOR	pDONR223
NCBI SYMBOL	VPS45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]