

C9orf9-chromosome 9 open reading frame 9 Gene

Reference: PTXBC012940

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf9; Mast; sperm acrosome-associated protein 9; Rsb66 homolog; sperm acrosome associated 9		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 9
GENE ID	11092
ORIGIN SPECIES	Human
ACCESSIONS	BC012940
VECTOR	pDONR223
NCBI SYMBOL	C9orf9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgaggtgaaagaatcccttcgagcatcgagcagaagtaagctcttcagcagcagctcacctcaccgccgtctggagcactgcagggagaagcccccagacaagatccggcccatctccagcattggacaggtgcagagctacatggaacactactgcaacagctccacagaccggcggttctgctcatgttcttgacatctgttcagagctgaataagctctgcagcactttgagccgtgcactctggcaccagtcaccaacaacctcctggagaaatgcaaaacctcgttagccaaagcaacgacttaagcagcctcagagcaaaataccctcatgatgtgtgtaaccacctcagctgtgacgagggccgggaaccactacggcggtgtgagcctcatccccctcatcctagactaatgaaagaatggatgcgccactccgagaagttgccgcaaggtgctgcagcagctgagtgagccccagggcgaccaggagagcaccaggggagccgctcctgccaggccataggacccagccagggccactaaacacaagtgtagacagctcacaagccagcctcaaacccaggggatgttcaaaaccctcgaggcctcctgtgggaaattgtaa