



Reference: PTXBC013007

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BTAF4; RAP74; TF2F1; TFIIF; general transcription factor IIF subunit 1; TFIIF-alpha; general transcription factor IIF 74 kDa subunit; general transcription factor IIF, polypeptide 1, 74kDa; transcription initiation factor IIF subunit alpha; transcription i		

For research use only.

GENE DESCRIPTION	general transcription factor IIF, polypeptide 1, 74kDa
GENE ID	2962
ORIGIN SPECIES	Human
ACCESSIONS	BC013007
VECTOR	pDONR223
NCBI SYMBOL	GTF2F1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]