

MRPS26-mitochondrial ribosomal protein S26 Gene

Reference: PTXBC013018

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf193; GI008; MRP-S13; MRP-S26; MRPS13; NY-BR-87; RPMS13; dJ534B8.3; 28S ribosomal protein S26, mitochondrial; 28S ribosomal protein S13, mitochondrial; S13mt; S26mt; serologically defined breast cancer antigen NY-BR-87; mitochondrial ribosomal protei		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S26
GENE ID	64949
ORIGIN SPECIES	Human
ACCESSIONS	BC013018
VECTOR	pDONR223
NCBI SYMBOL	MRPS26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctacgcgcgtgagccgctgggcgcgggacccgtgcaggccccgggcccctctggtgctgccagcgcgcggccgaagaccgccagaccgctggccaaatccaagatcgagcagtgaaacatgccgccgcgggtggaccctgcggagttctctgctgatggag
cggtaccagcactaccgccagaccgtgctgcgcctcaggatggagttctgttccgaggtgcagaggaagggtgcacagggccgagccgggggttctgctggagcgcgaagccctgaaggacgccgccgagcaccgcgagctgatggcctggaaccaggcggagaaccggcgc
tgcacgagctgcggatagcagagctgcgcgaggagcgggagcaggagcagcggcagggcgttggagcagcccgcaagccgaagaggtgcagccctgggcgcagcgaaggagcgggaagtgtcagctgcaggaagaggtgaaaaactcatcaccggagagaac
ctggagcagcgggtggaagcagcattggactcccggaagaactacaactggccatcaccagagagggcgtgtgtcaggccacaacgcagggaactcctag