



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

RILPL2-Rab interacting lysosomal protein-like 2 Gene

Reference: PTXBC013042

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RLP2; RILP-like protein 2; p40phox-binding protein; rab-interacting lysosomal protein-like 2; Rab interacting lysosomal protein like 2		

For research use only.

GENE DESCRIPTION	Rab interacting lysosomal protein-like 2
GENE ID	196383
ORIGIN SPECIES	Human
ACCESSIONS	BC013042
VECTOR	pDONR223
NCBI SYMBOL	RILPL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]