

FLJ14107-hypothetical protein FLJ14107 Gene

Reference: PTXBC013100

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BIN3 intronic transcript 1 (non-protein coding)		
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ14107
GENE ID	80094
ORIGIN SPECIES	Human
ACCESSIONS	BC013100
VECTOR	pDONR223
NCBI SYMBOL	FLJ14107
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcctccactggagaaggattttcatctgccaggtggcggggcccttagcctcgttaatcctgcttaacctgctagaacatactacaggagccactttccacgcagagcatccctgaaggaggggcttccagaatgctgaggattctgggtctcagactccgatcagatg
cactcagccctgggaltcagtgccccagccctgctgccagtttcttttcaactcaccatgagcataagttgggtgggaactcagcctgttcattcacttcatactgagttccactctgtgccaggcctgacctgtgctgccccaaagagcagggaacagcaaacctttctccac
cgtgggccttgagaaatgaatgcaggggggtggccacgggtag