

RHBDL2-rhomboid, veinlet-like 2 (Drosophila) Gene

Reference: PTXBC013103

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RRP2; rhomboid-related protein 2; rhomboid (veinlet, Drosophila)-like 2; rhomboid protease 2; rhomboid, veinlet-like 2; rhomboid-like protein 2; rhomboid like 2		

For research use only.

GENE DESCRIPTION	rhomboid, veinlet-like 2 (Drosophila)
GENE ID	54933
ORIGIN SPECIES	Human
ACCESSIONS	BC013103
VECTOR	pDONR223
NCBI SYMBOL	RHBDL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatctgaatatggggagagagatgaaagaagagctggaggaaaggagagaaatgagagagatgggggaggtaaaagatgggccaagagtaaaaaggtccacaggattgtctcaaatggatgctgcccgaagaggaacatactggagagagtaactgcttccgcctcccgtgtcatcatcctcagcctggccgagctggcagtggtttatttactatgctgtgtggaagcctcagaacagtgatcacgttgacacagcgcatttgagagtcctttatctacagtcctgagaagaggagggaagcctggaggtttatctcatcatgctggtacatgctgggttaa