



Reference: PTXBC013119

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPAL2; signal-induced proliferation-associated 1-like protein 2; SIPA1-like protein 2; SPA-1-like 2; signal induced proliferation associated 1 like 2		

For research use only.

GENE DESCRIPTION	signal-induced proliferation-associated 1 like 2
GENE ID	57568
ORIGIN SPECIES	Human
ACCESSIONS	BC013119
VECTOR	pDONR223
NCBI SYMBOL	SIPA1L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]